

Using cell free DNA (cfDNA) and RNA (cfRNA) in the diagnosis and monitoring of primary and metastatic CNS tumors - Abstract Code: INNV-08

Maher Albitar¹, Ahmad Charifa¹, Sally Agersborg¹, Omar Butt², Adam Albitar¹, Na Tosha N. Gatson³, and Santosh Kesari⁴.
¹Genomic Testing Cooperative; ²Washington University, St. Louis, MO; ³Indiana University, Indianapolis, IN; ⁴Saint John's Cancer Institute, Torrance, CA.



Background

In the current study, we expand our prior study and present results of sequential CSF samples from patients with suspected CNS neoplasm submitted for testing by LBx to rule out or confirm CNS tumor.

Methods

Each sample was centrifuged to isolate the cells. From the supernatant, cell-free DNA(cfDNA) and cell-free RNA (cfRNA) were extracted. From the cell pellet, cellular RNA (cRNA) was extracted. cfDNA was sequenced by next generation sequencing (NGS) using a targeted DNA panel of 302 genes. The cfRNA and cellular RNA were sequenced independently using a targeted RNA panel of more than 1,600 genes. From January 2024 to August 2025, 1397 consecutive samples were analyzed, 6 of which yielded inadequate cfDNA or cfRNA for analysis. Of the negative cases **35% were submitted with various diagnoses related to CNS including malignant neoplasm of unspecified brain or secondary neoplasm of central nervous system. 15% of negative cases were submitted with a diagnosis of lung cancer, 6% with a diagnosis of breast cancer and 4% with lymphoma.**

Results

Characteristics of the tested samples and the detected abnormalities.

Clinical Diagnosis	Number of samples	Number of detected SNV	Detected Abnormalities	Number	Percentage
ALL	14 (1.3%)	14	Total samples	1391	100.0%
AML	10 (0.9%)	22	Negative for SNV, CNV and clonality	231	16.6%
Brain	325 (30.2%)	1885	Positive for CNV	535	38.5%
Breast	138 (12.8%)	965	Pos for CNV only	14	1.0%
Carcinoma	62 (5.8%)	225	B-cell clonality	108	7.8%
Colorectal	3 (0.3%)	13	B-cell clonality only	2	0.1%
Endometrial/Ovarian	9 (0.8%)	72	Fusion genes	68	4.9%
Lung	138 (12.8%)	715	CAR-T cells	5	0.4%
Lymphoid	80 (7.4%)	447	Cases with Germlines	445	32.0%
Melanoma	6 (0.6%)	74	CHIP only	69	5.0%
Neuroendocrine	7 (0.7%)	33	17p(TP53) deletion	167	12.0%
Sarcoma	5 (0.5%)	25	TP53 Mutation+Deletion	96	6.9%
Other	278 (25.9%)	539	EGFR Amplification	30	2.2%
Total	1075	5029	ERBB2 amplification	31	2.2%


Detection of B-cell clonality in CNS lymphoma

Heavy chain	
CDR3	Counts
IGHV3-7	7197
IGHV3-62	54
IGHV7-51-2	2
IGHV3-63	2
IGHV3-65	1
IGHV7-65-1	1
Kappa	
CDR3	Counts
IGKV1-8	1040
IGKV2-26	143
IGKV1-12	1
IGKV7-3	1
IGKV1D-16	1
IGKV3-31	1

Tumor Load Vs Variant Allele Frequency

The bars represent the tumor load (number of mutant molecules in 1 mL of plasma). The VAF data suggests that the ESR1 mutation is in a subclone of that with the PIK3CA, CDKN2A and NF1. The tumor load is more representative and shows that one clone has PIK3CA, CDKN2A, NF1 and ESR1 mutations. The difference in VAF is due chromosomal copy number variation.


Gene	VAF %	Chr
PIK3CA	59.26	3
CDKN2A	48.12	9
NF1	41.95	17
PAK3	32.02	X
GRIN2A	29.81	16
ESR1	24.41	6
ARID2	22.93	12
KMT2C	12.9	7
AXIN2	11.02	17
BCOR	4.84	X




Most common mutations detected in various neoplasms

ALL	CEBPA	14.3%
AML	NPM1	20.0%
Brain	TP53	16.9%
Breast	TP53	37.0%
Carcinoma	TP53	22.6%
CRC	APC	33.3%
Endometrial/Ovarian	TP53	55.6%
Lung	TP53	29.7%
Lymphoid	KMT2C	18.8%
Melanoma	BRAF	50.0%
Neuroendocrine	KIT	28.6%
Sarcoma	CBL	40.0%
Other	TP53	22.7%

example of breast cancer with only chromosomal abnormalities and a second sample from a patient with EGFR gene amplification. Both had no detectable mutations



Detection of CAR-T cells in CSF



Conclusions

- Liquid biopsy of CSF is highly reliable in the diagnosis and monitoring of primary and Metastatic tumors.
- Clonality (B and T-cell) and chromosomal abnormalities should be a part of CSF evaluation
- RNA evaluation is required for detecting chromosomal translocations
- RNA can be used in evaluating ICAN Vs tumor in CSF.
- Quantifying the tumor load in CSF provides important information that are different from VAF.